

Duray[®] Intraoral Sensor Package

HOT DEAL

- Thinnest Sensor 4.4mm
- Durable and Reliable
Waterproof with IP68 protection and can be immersed in disinfectant liquids to prevent cross infections
- High Resolution Imaging
With a theoretical resolution of 25 lp/mm and through advanced CMOS technology, image quality is exceptional
- Easy To Use
Simple and easy-to-operate software with one-click smart installation and offers plug-and-play simplicity with high-speed USB 2.0 connection
- 24 Months Warranty



Duray Sensor Package [DS3001]
Reg. \$6,298 **Sale: \$5,250**

Size 1 [DS1001]
Reg. \$2,940 **Sale: \$2,450**

Size 2 [DS2001]
Reg. \$3,676 **Sale: \$3,125**



AirLight[®] M800 Plus

Generator + LED High Speed

- 23W of Torque
- 410,000 RPM
- Ceramic Ball Bearings
- Titanium Body
- 1 Year Warranty

Backend Code
4 Hole HP3008P
Beyes[®] HP6018P
W&H[®] HP3022P
Star[®] HP3054P
NSK[®] HP3043P
KAVO[®] HP3036P

Reg. \$548/ea **Net: \$411/ea**

BUY 3, GET 1 FREE!

Triple Spray
Direct-LED


Built-in Generator



PD Coupler
Non Optic, 4 Holes

Reg. \$152/ea **Net: \$114/ea**

BUY 3, GET 1 FREE!



Code: HP6008

Slimax-C Plus[®]

LED Curing Light System

- Powerful up to 2000 mW/cm2 (±10%)
- Built-in Digital Radiometer
- 3 Modes: Full Power, Ramping, Pulsing
- Rechargeable Lithium-ion Battery
- Cordless Design
- 5-Second Curing
- 1 Year Warranty

Color Code
Silver SM1003P-S
Blue SM1003P-B
Pink SM1003P-P

Reg. \$660/ea **Net: \$440/ea**



Maxso[®] M200E

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Single Spray
- Push Button
- Low Turbine Noise
- No Swivel
- 6 Months Warranty

Backend Code
4 Hole HP2105

Reg. \$128/ea **Net: \$110/ea**

BUY 5 FOR \$550!



Maxso[®] M200 Smart

Non-Optic High Speed Handpiece

- 19W of Torque
- 400,000 RPM
- Push Button Chuck
- Micro-tex Surface
- Triple Water Spray
- 1 Year Warranty

Backend Code
4 Hole HP2009

Reg. \$240/ea **Net: \$180/ea**

BUY 3, GET 1 FREE!



AirLight[®] M800X

Generator + LED High Speed

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Months Warranty

Backend Code
Beyes[®] HP6018X
Star[®] HP3054X
NSK[®] HP3043X
KAVO[®] HP3036X

Reg. \$1,275/ea **Net: \$637.50/ea**

BUY 1, GET 1 FREE!

Quattro Spray
Direct-LED

Built-in Generator



Maxso[®] Low Speed

E-Type Connection



1 Year Warranty **Mix and Match**

	Reg	Net
A. SL2026 1:1 Push Button (CA Burs)	\$176	\$132
B. SL2028 1:1 Latch Type (CA Burs)	\$147	\$110.25
C. SL2034 1:1 Screw-On	\$137	\$102.75
D. MT2010 Air Motor, 4 Hole, 20,000RPM	\$298	\$223.50
E. ST2010 1:1 Straight Nose Cone	\$114	\$85.50
F. SL1002 4:1 Nose Cone, For DPA	\$278	\$208.50

Mix & Match any 3. Promotional goods must be of equal or lesser value than the lowest priced item purchased.

AirLight® M800 Surgical 45

Non-Optic High Speed Handpiece

- Direct LED
- 22W of Torque
- 410,000 RPM
- 45° Head
- Rear air exhaust
- Triple water spray (Stream, not mist)
- Push Button
- 12 Months Warranty

Backend	Code
4 Hole	HP3048
Beyes®	HP6038
Star®	HP3055
NSK®	HP3044
KAVO®	HP3037

Reg. \$598/ea
Net: \$448.50/ea



BUY 3, GET 1 FREE!

Maxso® X200X

Non-Optic High Speed Handpiece

- 27W of Torque
- 430,000 RPM
- Instant Stop Capability
- Quattro Spray
- Titanium Body
- 25 Month Warranty

Backend	Code
NSK®	HP2066X
KAVO®	HP2060X

Reg. \$1,048/ea
Net: \$524/ea



BUY 2, GET 2 FREE!

EZ-Swivel®

Couplers



Beyes® PD 4 Hole
HP6008
Reg. \$152/ea
Net: \$101.33/ea



NSK® Type 4 Hole, Fiber-optic
HP3046
Reg. \$252/ea **Net: \$189/ea**

NSK® Type 6 Hole, Fiber-optic
HP2069
Reg. \$252/ea **Net: \$189/ea**



Kavo® Type 4 Hole,
Non-optic
HP2063
Reg. \$252/ea
Net: \$189/ea



Kavo® Type 6 Hole,
Fiber-optic
HP2064
Reg. \$252/ea
Net: \$189/ea

Maxso® X200

Fiber-optic High Speed Handpiece

- 23W of Torque
- Triple Spray
- Push Button Chuck
- Swivel
- 1 Year Warranty

Backend	Code
Kavo®	HP2060
NSK®	HP2066

Reg. \$480/ea **Net: \$360/ea**



BUY 3, GET 1 FREE!

Maxso® E800

Electric Micromotor System Package

- Powerful: 5.0 Ncm of torque
- Steady Performance
- Smooth Operation
- Convenient & Durable
- 3 Programmable Memory Settings
- New Touch Screen Interface
- Warranty: 25 months

Code: MT8001

Reg. \$1,726 **Sale: \$1,626**



Maxso® Low Speed

E-Type Connection, S20A-IS/E6

- Straight Nose Cone
- 1:1
- Internal Spray
- Fiber-Optic
- 1 Year Warranty

Code: ST2009

Reg. \$414/ea **Net: \$310.50/ea**



BUY 3, GET 1 FREE!

Maxso® X99L

1:5 Contra Angle

- Anti-Heat System
- Ceramic Ball Bearing
- ISO E-Type Connection
- Quattro Spray
- Fiber-optic
- 1 Year Warranty

Code: SL3025

Reg. \$992/ea **Net: \$744/ea**



BUY 3, GET 1 FREE!

Maxso® Low Speed

E-Type Connection, M29L-MH01

- Contra Angle
- Internal Spray
- Fiber-Optic
- Push Button
- CA Burs
- 1 Year Warranty

Code: SL3068

Reg. \$430/ea **Net: \$322.50/ea**



BUY 3, GET 1 FREE!