



Fiber-optic

Maxso®

X200 Surgical 45

It provides excellent accessibility, and visibility to make stress-free treatment a reality. X200 Surgical 45 features 25W of cutting power for faster preparation and rear exhaust to help increase patient safety. With its uniquely angled head, enhanced access is simple. Experience for yourself, the difference between ordinary and extraordinary.




Accessibility

45° Head design for perfect access to the molars.



Visibility

Stream water (not mist) for perfect bur cooling and visibility.



Hygiene

Rear air exhaust to prevent contamination into cutting area.



Power

Strong chuck holding force to prevent slippage and improve power with 25W torque.

Technical Data

Head Size (Width x Height)	10.8mm x 13.0mm
Head Angle	45°
Drive Air Pressure (PSI)	40-60
Power (WATT)	25
Speed (RPM)	400,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995
Ø / max. length (MM)	1.6 / 25
Chuck Type	Push Button
Water Port	Triple Water Jet

Order Info.

Kavo® Connection



Mini Head

Code **HP2062**

Model **X200-45/K**

NSK® Connection



Mini Head

Code **HP2068**

Model **X200-45/N**



Beyes
The Next Standard

Beyes Dental Inc.
55 Laurier Drive
Morrisburg, Ontario
K0C 1X0, Canada
Tel: 1-855-603-1888
Email: info@beyes.ca
Website: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands
Email: info@lotusnl.com

Your Dealer for Beyes Products



Scan the QR code
for more information

Printed in Canada
File #: ENB055R1
Rev.2/11.13.25

© Beyes Canada Inc 2025. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under the
Medical Device Single Audit Program (MDSAP).